

Supplementary Table S1 – AFLP primer combinations used for selective and pre-selective amplification.

Primer combinations – selective amplification:

Set A D2-EcoRI primer E-ACC / MseI primer M-AA

Set B D3-EcoRI primer E-ACT / MseI primer M-TT

Set C D4-EcoRI primer E-AGC / MseI primer M-AA

Pre-selective amplification primers:

EcoRI primer: 5'-GACTGCGTACCAATTCA-3'

MseI primer: 5'-GATGAGTCCTGAGTAAC-3'
