Supplementary Table S1 – AFLP primer combinations used for selective and pre-selective amplification.

Primer combinations – selective amplification:	
Set A	D2-EcoRI primer E-ACC / MseI primer M-AA
Set B	D3-EcoRI primer E-ACT / MseI primer M-TT
Set C	D4-EcoRI primer E-AGC / MseI primer M-AA
Pre-selective amplification primers:	
EcoRI primer: 5′-GACTGCGTACCAATTCA-3′	
MseI primer: 5′-GATGAGTCCTGAGTAAC-3′	